

Fluo-3, AM

Catalog Number	Packaging Size
C213	1 mg
C214	20×50 µg
C215	1 ml


Storage upon receipt: -20°C, protected from light

Introduction

Fluo-3, AM is a membrane permeant probe and can be passively loaded in cells by simple incubation. Once inside the cell, Fluo-3 becomes fluorescent in the presence of free Ca^{2+} . Fluo-3 AM is widely used in flow cytometry, confocal laser scanning microscopy and fluorescence spectroscopy. Fluo-3, AM has proved particularly valuable for determining GPCR function via cell-based high-throughput screening.

Fluo-3, AM is supplied as 1 mg solid (**C213**), special packing 20×50 µg (**C214**), and 1 mg/ml solution in DMSO (**C215**).

Specifications

Label:	Fluo-3	
Ex/Em:	488/525 nm	
Detection Method:	Fluorescent	
Molecular Formula:	$\text{C}_{51}\text{H}_{50}\text{Cl}_2\text{N}_2\text{O}_{23}$	
Molecular Weight:	1129.85	
CAS Number:	121714-22-5	
Storage Conditions:	-20°C, protect from light	
Shipping Condition:	Room Temperature	

Applications

Cell Loading Guideline

Note: The following protocol is provided as an introductory guide only. The detailed procedures can be found from literatures.^{1,2}

1. Prepare a Fluo-3 AM stock solution in anhydrous DMSO at 1-5 mM.
2. Dilute an aliquot of Fluo-3 AM stock solution (1-5 mM) to a final concentration of 1-5 µM in the buffered physiological medium of choice. Addition of the non-ionic detergent Pluronic R F-127 can assist in dispersion of the nonpolar Fluo-3 AM ester in aqueous media. This can be conveniently accomplished

by mixing the aliquot of Fluo-3 AM ester stock solution in DMSO with an equal volume of 20% (w/v) Pluronic in DMSO (Cat No. C021) before dilution into the loading medium, making the final Pluronic concentration about 0.02%.

3. The organic anion-transport inhibitors probenecid (1-2.5 mM) may be added to the cell medium to reduce leakage of the de-esterified indicator.
4. Cells are normally incubated with the Fluo-3 AM ester for 15–60minutes at 20–37°C. Exact loading concentration, time, and temperature will need to be determined empirically; in general it is desirable to use the minimum dye concentration required to yield fluorescence signals with adequate signal to noise. Subcellular compartmentalization, an inherent problem with the AM ester loading technique, is usually lessened by lowering the incubation temperature.
5. Before fluorescence measurements are commenced, cells should be washed in indicator-free medium (containing an anion transport inhibitor, if applicable) to remove any dye that is nonspecifically associated with the cell surface, and then incubated for a further 30 minutes to allow complete de-esterification of intracellular Fluo-3 AM ester.

High-Throughput Screening

Intracellular Ca²⁺ measurements in 96-well and 384-well microplates are an essential tool for high-throughput pharmacological screening. Cell samples in microplate wells are loaded with the AM ester form of the indicator using protocol basically similar to those described in ***Cell Loading Guideline***.

References:

1. Methods Cell Biol 40, 155 (1994);
2. *Cell Biology: A Laboratory Handbook, 2nd Edition*, J.E. Celis, Ed., Volume 3, pp 363–374, Academic Press (1998);