



NGS End Repair Module

Catalog Number: E112-1, E112-2

Table 1. Kit Components and Storage

Kit Component	E112-1 (20 rxn)	E112-2 (100 rxn)	Storage	Stability
End Repair Enzyme Mix	100 μ L	500 μ L	-20 °C, avoid repeated free-thaw	The product is stable for 12 months when stored as directed.
End Repair Reaction Buffer (10x)	200 μ L	1.0 mL		

Product Description

The NGS End Repair Module is used to convert 1–5 μ g of fragmented DNA to blunt-ended DNA having 5' phosphates, and 3'-hydroxyls. The module is optimized for use with the NGS dA-Tailing Module (Cat #E113-1, E113-2), and is part of the NGS DNA library prep workflow for Illumina sequencing.

Applications

- ❖ End repair of fragmented DNA.
- ❖ DNA sample preparation.

Protocol

- Assemble the following reaction in a microcentrifuge tube on ice:

Fragmented DNA	1-5 μ g
End Repair Reaction Buffer (10x)	10 μ L
End Repair Enzyme Mix	5 μ L
Nuclease-free water	to 100 μ L
Total volume	100 μ L
- Mix gently and spin down for a few seconds.
- Incubate at 20°C for 30 min.
- Purify DNA Sample using AMPure XP beads.