

## GST-Tag Mouse Monoclonal Antibody

Catalog Number: AT009

### Product Package

Cat. No.	Unit Size	Concentration	Storage
AT009	100 $\mu$ L	1 mg/mL	Short term: 4°C Long term: -20°C Avoid freeze/thaw cycle

### Introduction

GST-Tag antibody is part of the Tag series of antibodies, the best quality in the research.

**Storage Condition and Buffer:** 1 mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

**Application:** WB


**Reactivity:** All

**Dilution:** WB 1:1000~1:5000.

**Purification:** Affinity-chromatography.

**Specificity:** GST-Tag Mouse Monoclonal antibody detects endogenous levels of C-terminal, internal, and N-terminal GST-tagged proteins.

**Immunogen:** A synthesized peptide.



Western blot analysis of GST protein expression in (1) 293T cell lysate; (2) 293T cell transfected with GST protein lysate.

For Research Use Only. Not For Use In Diagnostic Procedures.  
[www.abpbiotech.com.cn](http://www.abpbiotech.com.cn)