


DMG-PEG2000

Catalog Number	Packaging Size
LP004-1	500 mg
LP004-2	1 g

Storage upon receipt: -20°C

DMG-PEG(2000) is a lipid excipient that has been used in combination with other lipids in the formation of lipid nanoparticles (LNPs). Formulations containing DMG-PEG(2000) have been used in the development of lipid nanoparticles for the delivery of mRNA-based vaccines.

For research use only.



Technical Information

Formal Name:	1,2-Dimyristoyl- <i>rac</i> -glycero-3-methoxypolyethylene glycol-2000
CAS Number:	160743-62-4
Purity:	>95%
Storage Condition:	-20°C
Shipping Condition:	Room Temperature